


No.P311

피크만에스피 HPLC Application

Coenzyme Q10

코엔자임 큐텐

1. Coenzyme Q10
(200 µg/mL)



[HPLC Conditions]

System	: NANOSPACE SI-2, OSAKA SODA, JAPAN
Column	: CAPCELL PAK C18 MG (4.6 ×150 mm, 5 µm), OSAKA SODA, JAPAN
Mobile phase	: ACN / THF / DW (55 : 40 : 5)
Flow rate	: 1.0 mL/min
Temperature	: 40°C
Detection	: PDA (280 nm)
Inj. vol.	: 10 µL
Sample dissolved in	: Hexane / EtOH (5 : 2)