


Flavonoid

CapcellPAK UG120 C18 4.6 x 250 mm, 5 μm



이동상 : A – DW / Acetonitrile / Phosphoric acid = 95 / 5 / 0.1

B – DW / Acetonitrile / Phosphoric acid = 5 / 95 / 0.1, Gradient

온도 : 40 °C

검출기 : PDA (254 nm)

Time (min)	Flow rate (μL/min)	B (%)	A (%)
Int.	1000	5	95
30	1000	25	75
35	1000	45	55
38	1000	45	55
40	1000	5	95
48	1000	5	95