



# Paeoniflorin

패오니플로린

CAPCELL PAK C18 KG

## 1. Paeoniflorin (50 µg/mL)



### [HPLC Conditions]

|              |                                                             |
|--------------|-------------------------------------------------------------|
| Column       | : CAPCELL PAK C18 KG (4.6 x 150mm, 3 µm), OSAKA SODA, JAPAN |
| Mobile phase | : 0.1% phosphoric acid in 14% ACN                           |
| Flow rate    | : 0.8 mL/min                                                |
| Temperature  | : 25°C, 35°C, 45°C, 55°C                                    |
| Detection    | : PDA (232nm)                                               |
| Inj. vol.    | : 10 µL (dissolved in 50% MeOH)                             |