



5-(3-Chlorophenyl)-3-hydroxypicolinic acid

5-(3-클로로페닐)-3-하이드록시피콜린산

CAPCELL PAK CN UG120

1. 5-(3-Chlorophenyl)-3-hydroxypicolinic acid
(100 µg/mL)



[HPLC Conditions]

Column	: CAPCELL PAK CN UG120 (4.6 × 250 mm, 5 µm), OSAKA SODA, JAPAN
Mobile phase	: A : 0.5% Phosphoric acid in DW B : 0.5% Phosphoric acid in ACN
Flow rate	: 1.0 mL/min
Temperature	: 30°C
Detection	: PDA (220 nm)
Inj. vol.	: 10 µL

